


Hill-Rom NP150 Viscoelastic

Self-adjusting, efficient and safe pressure ulcer prevention



Features and benefits

- Conforms to body shape and reduces interface pressure
- Welded cover with built-in carry handles
- Easy to clean and hygienic
- High density side sections designed for patient safety and improved egress
- Adapts to most bed frames and conforms to all bed positions
- X-ray translucent* and latex-free
- Optional integrated x-ray cassette sleeve assists in reducing patient discomfort when repositioning during portable x-ray procedures

* Models ASS095 and ASS099 only

Hill-Rom NP150 Viscoelastic

Temperature sensitive, viscoelastic foam mattress

The Hill-Rom NP150 Viscoelastic mattress has been specially developed to help prevent pressure ulcers. The temperature sensitive, viscoelastic foam moulds to the patient's shape, evenly distributes the weight and contours to the body in any position, increasing patient comfort and pressure redistribution. It is adaptable to most bed frames and conforms to all bed positions.

The fully enclosed mattress with its waterproof, yet vapour permeable, multi-stretch, bacteriostatic and fungostatic polyurethane cover is easy to clean and disinfect enhancing hygiene for patients and caregivers.



The Hill-Rom NP150 Viscoelastic offers 1 extended length

Standard size: 198 cm



Technical Specifications

Mattress (customised sizes available)

Length	198 or 203 cm
Width	85 or 90 cm
Height	14 cm

Warranty

Hill-Rom offers a 5 year warranty on the foam and a 2 year warranty on the cover. For details of the warranty provided please contact your local Hill-Rom representative.

Safe Working Load	250 kg
-------------------	--------

Recommended therapeutic weight	30 to 150 kg
--------------------------------	--------------

Seams	welded
-------	--------

Disinfection

Chlorine-based solutions can be used. Concentration must be less than or equal to 1,000 (or 10,000 ppm XT version).

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Hill-Rom NP 150 surface helps to prevent pressure ulcers in patients with no risk or up to moderate risk, weighing between 30 kg and 150 kg with a minimum height of 1.46 m, in a range of healthcare environments.

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN116303-04, 20 April 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom